

ClearSkin Pro 7x7 Pixel Tip

Indications:

- Aged/Photodamaged Skin
- Scars, including Acne Scars

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small section of the intended treatment area with the lowest recommended fluence before the first treatment session. Assess skin condition after 15-30 min.
- All personnel within the laser operating area must wear appropriate eye protection.
- It is recommended to use topical anesthetic cream on the treatment area before treatment.

Treatment:

- Position the applicator perpendicular to skin, cover treatment area with minimal overlap.
- Endpoint: moderate erythema.
- Number of treatments: 4-6.
- Treatment intervals: 2-3 weeks.



Recommended Treatment Parameters – Aged and Photodamaged Skin/Scars:

Fitzpatrick Skin Type	Energy (mJ/P)	Pulse Repetition Rate (Hz)	# Stacks	# Passes
I – III	2000 – 3000	1 – 2	1 – 2	2 – 3
IV	1500 – 2500	1 – 2	1 – 2	1 – 2
V-VI	1500 – 2000	1 – 3	1	1

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema apply Trolamine, Cicaplast or equivalent skin calming products.

ClearSkin PRO CVD Tip (2 – 10 mm)

Indications:

- Acne Vulgaris
- Scars, including Acne Scars

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area of the intended treatment area with the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.

Treatment:

- Position the applicator perpendicular and treat lesion by lesion. Detach vacuum by tilting applicator sideways (do not pull away).
- Number of treatments: 4 – 6.
- Treatment intervals: 2 – 3 weeks.



Recommended Treatment Parameters – Acne Vulgaris/Scars:

Fitzpatrick Skin Type	Spot Size (mm)	Energy (mJ/P)	Pulse Repetition Rate (Hz)	# Stacks	# Passes
I – III	3	300	1	1 – 2	1
	4 – 5	300 – 500			
	6-7	500- 800			
	8 – 10	800-1000			
IV – VI	4 – 5	300	1	1 – 2	1 – 2
	6-7	400-700			
	8 – 10	600-800			

Vacuum – On

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30)
- If treatment area shows erythema apply Trolamine, Cicaplast or equivalent skin calming products.