

Alma Iris Dye-VL SVL (In Motion) Mode

Indications:

- Skin Vascular Lesions
- Skin Pigmented Lesions
- Aged/Photodamaged Skin

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on the intended treatment area using the lowest recommended fluence before the first treatment session. Assess skin condition after 30 min for skin types I – III and after 24 – 48 hours for skin type IV – V.
- All personnel within the AFT operating area must wear appropriate eye protection.

Treatment:

- Apply ultrasound gel – coupling medium.
- Position the applicator perpendicular to skin.
- Endpoint: mild erythema, darkening of pigmented lesions, change in vessel coloration (reddening, darkening, or bleaching)
- Number of treatments: 3 – 6 sessions.
- Treatment intervals: 2 – 3 weeks.



Recommended Treatment Parameters – Pigmented and Vascular Lesions, Photodamaged Skin: SVL Mode (In-Motion):

Fitzpatrick Skin Type	Cycle Duration Time (s)	Fluence (J/cm ²)	Total Energy (kJ)*	Cooling (%)
I – II	30	3 – 4	3.2 – 4.3	100
III	30	2 – 3	2.1 – 3.2	100
IV	30	2	2.1	100

* Energy is calculated for 10x10 cm (100 cm²) area.

**Recommended Treatment Parameters – Pigmented and Vascular Lesions, Photodamaged Skin: SVL Mode (In-Motion):
Small Treatment Area (3-25 cm²):**

Fitzpatrick Skin Type	Cycle Duration Time (s)	Fluence (J/cm ²)	Total Energy (kJ)*	Cooling (%)
I – II	3	3 – 4	0.3 – 0.4	100
III	3	2 – 3	0.2 – 0.3	100
IV	3	2	0.2	100

* Energy is calculated for up to 5x5 cm (25 cm²) area.

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.

Alma Iris Dye-VL VL (Stationary) Mode

Indications:

- Skin Vascular Lesions
- Skin Pigmented Lesions
- Aged/Photodamaged Skin

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on the intended treatment area using the lowest recommended fluence before the first treatment session. Assess skin condition after 30 min for skin types I – III and after 24 – 48 hours for skin type IV .
- All personnel within the operating area must wear appropriate eye protection.

Treatment:

- Apply ultrasound gel – coupling medium.
- Position the applicator perpendicular to skin.
- Pigmented lesions – cover treatment area with minimal overlap.
- Endpoint: mild erythema, darkening of pigmented lesions, change in vessel coloration (reddening, darkening, or bleaching)
- Number of treatments: 3 – 6 sessions.
- Treatment intervals: 2 – 3 weeks.



Recommended Treatment Parameters – Pigmented and Vascular Lesions, Photodamaged Skin: VL Mode (Stationary)

Fitzpatrick Skin Type	Pulse Width (ms)	Fluence (J/cm²)	Cooling (%)
I – III	10, 12	8 – 12	75/100
IV	12, 15	6 – 8	75/100

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.