

Alma SupErb 7x7 Pixel Tip

Indications:

- Aged/Photodamaged Skin
- Skin Pigmented Lesions
- Scars, including acne scars

Pre-Treatment:

- Recommended prophylactic oral antiviral herpes simplex Virus 1 (HSV-1), 24 hours prior to laser treatment.
- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area in the intended treatment area before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.



Treatment:

- Position the applicator perpendicular to skin, cover treatment area with minimal overlap.
- Can be used with the Roller endpiece.
- Use smoke evacuator and facial mask.
- Endpoint: epidermal ablation in a 7x7 pixel pattern.
- Number of treatments: 2 – 4 sessions.
- Treatment intervals: 3 – 4 weeks.

Recommended Treatment Parameters – 7x7 Pixel Tip

Fitzpatrick Skin Type	Pulse width (ms)	Energy* (mJ/P)	Frequency (Hz)	# Stacked Pulses**	# Passes
I – III	1.0-1.5	1600 – 2400	2	2 – 3	1 – 2
IV	0.75	1000-1600	2	2 – 3	1 – 2
V	0.25-0.5	600-1000	2	2 – 3	1 – 2

* Once the trigger is pressed, a one second preheating phase occurs before the first pulse.

** Treatment Intensity is defined by the number of stacked pulses.

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30)
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products

Alma SupErb 7x1 Pixel Tip

Indications:

- Aged/Photodamaged Skin
- Skin Pigmented Lesions
- Scars

Pre-Treatment:

- Recommended prophylactic oral antiviral herpes simplex Virus 1 (HSV-1), 24 hours prior to laser treatment.
- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area of the intended treatment area using the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.



Treatment:

- Position the applicator perpendicular to skin, cover treatment area with minimal overlap.
- Can be used with the Roller endpiece.
- Use smoke evacuator and facial mask.
- Endpoint: epidermal pixel ablation.
- Number of treatments: 2 – 4 sessions.
- Treatment intervals: 3 – 4 weeks.

Recommended Treatment Parameters – 7x1 Pixel Tip

Fitzpatrick Skin Type	Treatment Intensity	Pulse width (ms)	Energy (mJ/P)	# Passes *
I – III	Mild	0.25 – 0.75	<600	1 – 2
IV-V	Mild	0.25 – 0.75	<600	1
I – III	Moderate	0.5 – 1.0	<800	1 – 2
IV – V	Moderate	0.5 – 1.0	<800	1
I – III	Aggressive	0.75 – 1.5	<1200	1 – 2
IV -V	Aggressive	0.75 – 1.5	<1200	1

*- 1 Pass refers to treatment over area in a crosshatched pattern.

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.

Alma SupErb CVD Tip (1 – 6 mm)

Indications:

- Skin Pigmented Lesions
- Warts

Pre-Treatment:

- Recommended prophylactic oral antiviral herpes simplex Virus 1 (HSV-1), 24 hours prior to laser treatment.
- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area in the intended treatment area using the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.



Treatment:

- Position the applicator perpendicular to skin, cover treatment area with minimal overlap.
- Use smoke evacuator and facial mask.
- Endpoint: epidermal ablation/detachment of skin tags/ lesion resurfacing.
- Number of treatments: 1 – 3 sessions.
- Treatment intervals: 4 – 6 weeks.

Recommended Treatment Parameters – 1 – 6 mm CVD Tip:

Fitzpatrick Skin Type	Spot Size (mm)	Pulse Width (ms)	Energy* (mJ/P)
I – V	1 – 2	0.25 – 0.75	<1000
	3 – 4	0.5 – 1.0	<2000
	5 – 6	1.0 – 1.5	<2500

* – Once the trigger is pressed, a one second preheating phase occurs before the first pulse.

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products or dress the wound if necessary.

Alma SupErb CVD Tip (8 mm), Thermal Mode

Indications:

- Aged/Photodamaged Skin

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on the intended treatment area using the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.



Treatment:

- Position the applicator perpendicular to skin, move quickly within the grid area until total energy is reached.
- Endpoint: erythema.
- Number of treatments: 2 – 4 sessions.
- Treatment intervals: 3 – 4 weeks.

**Recommended Treatment Parameters – 8mm CVD Tip
(Thermal Mode):**

Fitzpatrick Skin Type	Treatment Area	Energy (mJ/P)	Total Energy (J)
I – V	Perioral	200 – 250	500 – 600
	Cheek (one side)	200 – 250	800 – 1000
	Crows feet	100 – 150	300 – 350
	Infraorbital	100 – 150	200 – 250

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.