

ClearLift Pro 5x5 Pixel Tip

Indications:

- Aged/Photodamaged Skin

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area of the intended treatment area with the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.

Treatment:

- Position the applicator perpendicular to skin, move slowly within the grid area until total energy is reached.
- Can be used with the Roller endpiece.
- Endpoint: Mild erythema.
- Number of treatments: 3 – 5 sessions.
- Treatment intervals: 2 – 3 weeks.



5x5 Pixel Tip Recommended Treatment Parameters – Aged/ Photodamaged Skin:

Fitzpatrick Skin Type	Indication	Focal Depth	Pulse Repetition Rate (Hz)	Energy (mJ/Pulse)	Total Energy (J)*
I-V	Aged Skin	(-2),(-1), 0	3-4	1800-2400	1000
V-VI			4-6	1200-1800	1000
I-IV	Photodamaged Skin / pigmentation	(+1), (+2)	4-6	1200-1800	600-800

* The total energy is recommended for 5x4 cm grids (20 cm²)

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30)
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.

ClearLift Pro CVD QS 1064 Tip

Indications:

- Unwanted Tattoo – Blue, Black, and Grey pigments
- Skin Pigmented lesions

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area in the intended treatment area with the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.



Treatment:

- Position the applicator perpendicular to the skin, cover tattoo/lesion with minimal overlaps.
- Endpoint: Frosting, skin elevation.
- Perform 1 – 2 test spots at lowest fluence and evaluate the skin's response.
- Number of treatments: 8 – 12 sessions for tattoo removal, and 3 – 6 sessions for pigmented lesions.
- Treatment intervals: 8 – 12 weeks for tattoo removal, and 3 – 4 weeks for pigmented lesions.

Recommended Treatment Parameters – Unwanted Tattoo and Skin Pigmented Lesions:

Fitzpatrick Skin Type	Spot Size (mm)	Energy (mJ/ Pulse)	Pulse Repetition Rate (Hz)
I – IV	3	400	1-2
	4	400 – 800	1-2
	5	400-1200	1-2
	6-7	400-1600	1-2
	8-9	800-1600	1-2
	10	1200-1600	1-2
V-VI	4-5	400	1-2
	6	400-1200	1-2
	7	400-1200	1-2
	8	800-1200	1-2
	9-10	800-1200	1-2

Adjust the energy according to the session number and patient's tolerance.

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.

ClearLift Pro CVD QS 532 Tip

Indications:

- Unwanted Tattoo – red, orange, and yellow tattoo pigments
- Superficial Skin Pigmented lesions

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on a small area of the intended treatment area at **lowest fluence** before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the laser operating area must wear appropriate eye protection.

Treatment:

- Position the applicator perpendicular to the skin, cover tattoo/lesion with minimal overlaps.
- Endpoint: Frosting, skin elevation
- Perform 1 – 2 test spots and evaluate the skin's response.
- Number of treatments: 8 – 12 sessions for tattoo removal, and 1 – 2 sessions for pigmented lesions.
- Treatment intervals: 8 – 12 weeks for tattoo removal, and 3 – 4 weeks for pigmented lesions.



Recommended Treatment Parameters – Unwanted Tattoo and Skin Pigmented Lesions (Superficial):

Fitzpatrick Skin Type	Spot Size (mm)	Energy (mJ/ Pulse)	Pulse Repetition Rate (Hz)
I – III	3 – 5	400	1-2
	6-10	400-800	1-2
IV	5 - 10	400	1-2

Adjust the energy according to the session number and patient's tolerance.

Post – Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30)
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.