

Alma Iris VL/PL

Indications:

- Skin Vascular Lesions

Pre-Treatment:

- Clean treatment area; remove all make-up or moisturizing cream.
- Always perform a skin test on the intended treatment area using the lowest recommended fluence before the first treatment session. Assess skin condition after 15 – 30 min.
- All personnel within the operating area must wear appropriate eye protection.

Treatment:

- Apply ultrasound gel – coupling medium.
- Position the applicator perpendicular to skin, cover treatment area with minimal overlap in a stationary technique.
- Endpoint: mild erythema, darkening of pigmented lesions, change in vessel coloration (reddening, darkening, or bleaching)
- Number of treatments: 3 – 6 sessions.
- Treatment intervals: 2 – 3 weeks.



Recommended Treatment Parameters – Vascular Lesions:

Fitzpatrick Skin Type	Pulse Width (ms)	Fluence (J/cm ²)	Cooling (%)
I-III	10, 12	14 – 23	75/100
IV	12, 15	12 – 15	75/100
V	15	8 – 10	75/100

Post-Treatment:

- Avoid sun exposure after and in between treatments, use UVA/UVB sunscreen (SPF>30).
- If treatment area shows erythema, apply Trolamine, Cicaplast or equivalent skin calming products.